

ULTRAcel

Multi-layered non surgical skin tightening

ULTRAcel works on all layers of your skin from the inside out. The combination of treatment options used will provide the ultimate results for your skin



UNTREATED

TREATED - IMMEDIATE EFFECT

ULTRAcel

Multi-layered non surgical skin tightening

New & Exclusive

ULTRAcel is a revolutionary, no downtime, triple facelifting technology, which uses the powerful but gentle ultrasound, focused radio frequency and fractionalized microneedling. These three techniques combine to lift not just the skin, but the deep muscle layer underneath, which was previously only reached by surgical facelifts.

Together they create the most complete non-surgical facelifting technology:

The Ultra-lift



For more information on ULTRAcel:
www.ultraceluk.com
+44 (0) 1481 736837

ULTRAcel

"The most advanced multi-layered skin tightening technology."

Amazing results after treatment"



*SPLIT FACE IMAGE TAKEN AFTER ONE TREATMENT.
Photography and quotes courtesy of Healthxchange Pharmacy Ltd & Dr Rita Rakus


www.ultraceluk.com

What is ULTRAcel?

ULTRAcel delivers 3 different modalities through one very clever system; Grid Fractional Radio frequency (GFR), Fractional Radio frequency Microneedling (FRM) and High Intensity Focused Ultrasound (FUS).

These modalities are used in combination or as stand-alone treatments. ULTRAcel achieves a combination of short and long term results with minimal or no downtime.

Depending on the treatment chosen by you and your practitioner, you can resume your normal daily routine the same day.



How does it work?

The radio frequency delivers heat, which affects collagen in the deeper layers of the skin. This causes deep structures of the skin to tighten, and over time, new and remodelled collagen is produced.

Ultrasound waves are then beamed to a controlled targeted depth of the skin where thermal heat energy tightens and contracts underlying tissue, right down to the muscle layer.

Thermal deep tissue treatments stimulate the growth of new healthy collagen for prolonged and more intensive tightening.

What will ULTRAcel do for me?

ULTRAcel is the ultimate treatment for working on and within your skin producing exceptional results with little or no downtime. ULTRAcel can help to:

- Reduce sagging by tightening your skin, typically jowls & neck
- Reduce fine lines & wrinkles
- Improve overall skin texture & tone
- Restore a youthful & smoother appearance
- Diminish large pores
- Treat acne & acne scarring

Additionally it works on a layer of the tissue that is generally only improved by surgery, the SMAS (Superficial Muscular Aponeurotic System).



“ I am absolutely thrilled with the results of my ULTRAcel treatment which I received yesterday. I can honestly say that I was rendered completely speechless as the result was as you said – immediate. ”

Mrs N, from Surrey

Frequently asked questions:

When will I see results?

The results from these treatments are variable, depending on the original quality of the skin and degree of damage to be treated. The skin will begin to feel smoother and tighter even after the first treatment, results from the treatment will continue to improve for up to six months, increasing with additional treatments.

How many treatments will I need?

Treatments are specific to you and your Medical Professional will discuss what it is you would like to achieve and how many treatments it will take for you to reach your goal. Most patients generally take between 1-3 treatments.

What aftercare routine should I follow?

It is recommended to use a medically proven skin care regime such as Obagi. Always use a sunscreen (a minimum of SPF30) to avoid sun damage.

How long before I can return to normal activities?

With GFR and FUS you can return to your normal activities immediately with minimal to no downtime. The FRM treatment may take between 3 hours and 3 days before normal activities are resumed.

What side effects can I expect?

Side effects are rare, have a chat with your doctor if you have any concerns.

